

 **UTM**<sup>®</sup> UNIVERSAL TRANSPORT MEDIA

COLLECTION, TRANSPORT AND PRESERVATION OF CLINICAL SPECIMENS FOR VIRAL MOLECULAR DIAGNOSTIC TESTING



ROOM TEMPERATURE  
STABLE



CLEARED FOR  
USE WITH VIRAL  
MOLECULAR ASSAYS




MAINTAINS VIABILITY  
OF ORGANISMS FOR  
UP TO 48 HOURS




CONTAINS  
ANTIBIOTICS TO  
INHIBIT BACTERIAL  
AND FUNGAL FLORA  
IN PATIENT SAMPLES

# UTM™ Collection Quick Guide



## COLLECT

Collect the patient sample using the swab. To limit the risk of contamination, be sure that the swab tip only touches the collection site and avoid touching the swab below the breakpoint.



## SNAP

Remove cap from the tube and insert the swab to the bottom of the tube.

Holding the tube away from face, grasp the end of swab shaft and bend at a 180-degree angle to break at the breakpoint. If needed, gently rotate the swab shaft to complete the breakage.

Screw the cap on tightly to prevent leakage.









## SEND

Identify tube with patient information and send to laboratory.

*Specimen collection should be performed by health care personnel who have completed training and demonstrated competency. Always read the manufacturer's package insert for specific instruction regarding specimen collection and transport for the type of test kit being used.*




## Product Description

	DESCRIPTION	CAT. NO.	PACKAGING	
KITS	<b>Minitip Size Flocked Swab</b> with Breakpoint Packaged with 3 mL UTM® Medium in 16x100mm Skirted Tube with Plastic Red Cap		3C057N	50/Pack; 6 Packs/Case
	<b>Minitip Size Foam Swab</b> with Breakpoint Packaged with 3 mL UTM® Medium in 16x100mm Skirted Tube with Plastic Red Cap		3C058N	50/Pack; 6 Packs/Case
	<b>Regular Size Polyester Swab</b> with Breakpoint Packaged with 3 mL UTM® Medium in 16x100mm Skirted Tube with Plastic Red Cap		3C059N	50/Pack; 6 Packs/Case
TRANSPORT MEDIA	3 mL UTM® Medium in 16x100mm Tube - Bulk		3C047N	50/Pack; 6 Packs/Case
SWABS	Single <b>Minitip Size Foam Swab</b> with Breakpoint		1C061S01	100/Pack; 10 Packs/Case
	Single <b>Regular Size Polyester Swab</b> with Breakpoint		1U054S01	100/Pack; 10 Packs/Case



 Made in the USA

For more resources visit [copanusa.com/sample-collection-transport-processing/utm-viral-transport/](http://copanusa.com/sample-collection-transport-processing/utm-viral-transport/).

 COPANUSA.com  
 800.216.4016  
 @COPANUSA

 info@copanusa.net  
 facebook.com/COPAN.USA