Abbott

CLEARVIEW™
PBP2a SA
CULTURE COLONY TEST

CLEARVIEW™ PBP2a SA CULTURE COLONY TEST

For *in vitro* Use Only

Rx Only

INTENDED USE

The Clearview™ PBP2a SA Culture Colony Test is a qualitative, *in vitro* immunochromatographic assay for the rapid detection of penicillin-binding protein 2a (PBP2a) in isolates identified as *Staphylococcus aureus* as an aid in identifying methicillin-resistant *Staphylococcus aureus* (MRSA).

SUMMARY and EXPLANATION of TEST

Early detection of methicillin-resistant Staphylococcal infections is known to reduce health care costs, critical in efforts to decrease patient morbidity and mortality, reduce empiric use of vancomycin and permits cost-effective decisions for optimal patient management.¹

Conventional methods of antimicrobial susceptibility testing can take over 24 hours, and *mecA* gene identification is expensive and time-consuming to perform. However, penicillin-binding protein 2a (PBP2a) detection has the advantage over *mecA* in identifying strains that not only harbor the *mecA* gene but also produce the protein that confers resistance to methicillin.

The Clearview PBP2a SA Culture Colony Test provides a simple and rapid method to detect the PBP2a protein found in isolates of methicillin-resistant *Staphylococcus aureus*.

PRINCIPLES of PROCEDURE

The Clearview PBP2a SA Culture Colony Test is a rapid immunochromatographic membrane assay that uses highly sensitive recombinant monoclonal antibody fragments (rFabs) to detect the PBP2a protein directly from bacterial isolates. The rFab and a control protein are immobilized onto a nitrocellulose membrane as two distinct lines and combined with a sample pad, a pink/purple conjugate pad, and an absorption pad to form a test strip.

Isolates are sampled directly from the culture plate and eluted into an assay tube containing Reagent 1. Reagent 2 is then added and the test strip is placed in the assay tube. Results are read visually at 5 minutes.

REAGENTS and MATERIALS

Materials Provided in the CLEARVIEW™ PBP2a SA Culture Colony Test Kit

- **Test Strips:** A test strip covered with a plastic overlay with directional arrows to indicate testing position.
- **Reagent 1:** A clear, blue alkaline solution. 
- **Reagent 2:** A clear, slightly acidic solution containing sodium azide buffer and surfactants. 
- **Assay Tubes**
- **Test Racks**

MATERIALS REQUIRED, but NOT PROVIDED:

Clock, timer or stopwatch, bacteriological loops, vortex mixer and external positive and negative control strains.

EXTERNAL POSITIVE and NEGATIVE CONTROLS

The recommended controls are listed below. Alternatively, laboratories may use resistant and sensitive *S. aureus* strains routinely used, providing control results are as expected.

Positive control: *Staphylococcus aureus*, ATCC # 43300

Negative control: *Staphylococcus aureus*, ATCC # 25923

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. If refrigerated, allow all kit components to equilibrate to room temperature (15-30°C) before use.
3. Leave test strip sealed in its foil pouch until just before use.
4. Avoid skin and eye contact with reagents and test strip.
5. Do not use kit past its expiration date.
6. Do not interchange or mix components from different kit lots.
7. As with all bacterial samples, controls and test strips may contain pathogenic organisms; handle with appropriate precautions and dispose of materials safely in biohazard waste receptacles.
8. The Clearview PBP2a SA Culture Colony Test should be performed only on isolates of *Staphylococcus aureus*.
9. Reagent 1 contains sodium hydroxide
DANGER: Causes severe skin burns and eye damage. 
10. Reagent 2 contains sodium azide and surfactants.
WARNING: Causes serious eye irritation. 
11. Safety Data Sheets for this product are available upon request.
12. Follow your national, regional, and local ordinances accordingly for waste disposal regulations.

STORAGE and STABILITY

Store kit components at room temperature or under refrigeration (2-30°C).

The Clearview PBP2a SA Culture Colony Test Kit and reagents are stable until the expiration dates marked on their outer packaging.

QUALITY CONTROL

Daily Quality Control:

Clearview PBP2a SA Culture Colony Test has built-in positive and negative procedural controls. For quality control, Abbott suggests that you record these controls for each test run.

Procedural Controls:

- A. The appearance of a pink/purple line at the “control line” position can be considered an internal positive procedural control. If capillary flow has occurred, this line will always appear.
- B. In comparison to the color of the control line, the background color on the test strip should be white within 5 minutes.

EXTERNAL POSITIVE and NEGATIVE CONTROLS:

Good laboratory practice suggests the use of positive and negative controls to ensure that:

- test reagents are working, and
- the test is correctly performed.

For each new lot, external positive and negative controls should be tested. These will monitor the entire assay. The results of these control samples should be recorded. Please refer to the Sample Preparation Procedure section for further instructions.

The recommended controls are listed below. Alternatively, laboratories may use resistant and sensitive *S. aureus* strains routinely used, providing control results are as expected.

Recommended Control Strains:

Positive control: *Staphylococcus aureus* ATCC # 43300

Negative control: *Staphylococcus aureus* ATCC # 25923

Other controls may be tested in order to conform with,

- local, state, and/or federal regulations,
- accrediting groups, and/or,
- your laboratory's standard Quality Control procedures.

If the correct control results are not obtained, do not report results. Contact Technical Support (please refer to the Ordering and Contact Information section).

SPECIMEN COLLECTION and HANDLING

Specimens are bacterial isolates of *Staphylococcus aureus*. The use of fresh (<24 hours) cultures is recommended. The performance of the Clearview PBP2a SA Culture Colony Test has not been established for use with refrigerated specimens.

Culture Media

S. aureus colonies may be tested from any of the following culture media:

- Tryptone Soy Agar (Tryptic Soy Agar) with 5% sheep blood (TSA blood)
- Columbia Agar with 5% sheep blood
- Mueller Hinton Agar

Control Procedures

1. Subculture the control strain onto a culture plate. Incubate the plate overnight at 33-35°C for 18-24 hrs.
2. Follow Assay Procedure below.

Assay Procedure

If refrigerated, allow reagents and test strips to equilibrate to room temperature (15-30°C) before testing.

The test can be performed from well-isolated colonies on the primary plate if there is sufficient growth, or from a subculture of the isolate.

1. Holding the dropper bottle vertically, add two drops of Reagent 1 to an assay tube.
2. Take one heaped 1 µl bacteriological loop (a heavy inoculum) of sample from well-isolated colonies on the culture plate, place into the tube and thoroughly mix.
3. Holding the dropper bottle vertically, add two drops of Reagent 2 to the tube.
4. Vortex briefly. The blue solution must turn a clear color (if the color does not change, add one more drop of Reagent 2 and mix until the sample turns clear).
5. Insert the test strip into the assay tube with arrows pointed downward.
6. At five (5) minutes, withdraw the test strip from the tube and read the assay result.


RESULT INTERPRETATION

For help in reading the test strip, refer to the Result Interpretation Graphic.

For a **Negative Sample**, a PINK/PURPLE Control Line appears in the top half of the test strip. No other line appears.

For a **Positive Sample**, the PINK/PURPLE Control Line appears AND a second PINK/PURPLE Sample Line appears below it in the bottom half of the test strip. Any Sample Line, even when very faint, is positive.

A test is **Invalid** if the PINK/PURPLE Control Line does not appear, whether a Sample Line is present or not. Repeat invalid tests with a new test strip. Call Abbott Technical Support if the problem persists.

Reporting of Results

Result	Suggested Report
Positive	Positive for PBP2a
Negative	Negative for PBP2a
Invalid	Indeterminate; do not report. Re-test sample and contact Abbott Technical Service Support if problem persists.

CLINICAL PERFORMANCE

The clinical performance of the Clearview PBP2a SA Culture Colony Test was established in a multi-center clinical study conducted in 2013 at three (3) geographically diverse laboratories.

A total of 454 *S. aureus* samples were evaluated in the Clearview PBP2a SA Culture Colony Test, compared to results of 30 µg cefoxitin disk diffusion and interpreted according to CLSI standards. Clearview PBP2a SA Culture Colony Test performance versus cefoxitin disk diffusion, including 95% confidence intervals and stratified by plate type, is provided in Table 1.

All positive and negative daily controls generated the expected results.

Table 1: CLEARVIEW™ PBP2a SA Culture Colony Test Performance vs. Cefoxitin (30 µg) Disk Diffusion in *S. aureus* Isolates: Results by Plate Type

Plate Type	Sensitivity	95% C.I.	Specificity	95% C.I.
Primary Plate ¹	100.0% (129/129)	(97.1, 100.0)	98.5% (134/136)	(94.8, 99.6)
Tryptic Soy Agar with 5% sheep blood	99.1% (213/215)	(96.7, 99.8)	99.2% (237/239)	(97.0, 99.8)
Columbia Agar with 5% sheep blood	98.6% (212/215)	(96.0, 99.5)	99.2% (237/239)	(97.0, 99.8)
Mueller Hinton with 30 µg cefoxitin induction	99.1% (213/215)	(96.7, 99.8)	99.6% (238/239)	(97.7, 99.9)

1: The Clearview PBP2a SA Culture Colony Test was performed from primary plates at 2 out of 3 clinical sites. Primary plates were either Tryptic Soy Agar or Columbia Agar, with the exception of two samples of unknown plate type.

ANALYTICAL PERFORMANCE

Analytical Reactivity and Specificity

162 strains of methicillin-resistant *Staphylococcus aureus* (MRSA) and 112 strains of methicillin-sensitive *Staphylococcus aureus* (MSSA) were tested with Clearview PBP2a SA Culture Colony Test with expected results. These bacterial strains were obtained from the Network on Antimicrobial Resistance in *Staphylococcus aureus* (NARSA), American Type Culture Collection (ATCC) and a collection of strains from the Department of Infectious Disease Epidemiology of the Imperial College in London, England. Please note that isolates tested were not cultured using Columbia agar plates. All positive and negative daily controls generated the expected results.

REPRODUCIBILITY STUDY

A study of Clearview PBP2a SA Culture Colony Test was conducted at 3 separate sites using panels of blind coded specimens containing negative and positive samples in duplicate for each panel. Participants each tested one panel on 5 different days. There was 100.0% (600/600) agreement with expected test results. There were no significant differences within run (replicates tested by one operator), between run (5 different days), between sites (3 sites), or between operators (6 operators). Please note that isolates tested were not cultured using Columbia agar plates. All positive and negative daily controls generated the expected results.

SYMBOLS

 Consult instructions for use	 Authorized Representative in the European Community	 <i>In vitro</i> diagnostic medical device
 Catalog number	 CE Mark	 Manufacturer
 Temperature limitation	 Hazard Pictogram. See precautions.	 Kit contains sufficient materials for 25 tests
 Do not reuse	 Batch code	 Use by
Rx Only Prescription Only (Applies to US only)		

ORDERING and CONTACT INFORMATION

893-000 Clearview PBP2a SA Culture Colony Test

US: 1 877 441 7440

OUS: +1 321 441 7200

TECHNICAL SUPPORT

Advice Line

Further information can be obtained from your distributor, or by contacting Abbott Technical Support on:

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© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
Abbott Diagnostics Scarborough, Inc.