



Instructions for Use

COMPACT DRY™ BC

Cat. no. XBC100

Compact Dry™ BC, 60x75mm Tray with 10x60mm Well

100 trays/box

INTENDED USE

Hardy Diagnostics Compact Dry™ BC is a rapid, ready-to-use chromogenic test method recommended for the presumptive detection of *Bacillus cereus* in raw materials, finished products, or on environmental surfaces pertaining to food and other related industries.

This product is not to be used for the diagnosis of human disease.

SUMMARY

Bacillus cereus is a gram-positive, spore-forming microorganism common to soil, water, or other environmental samples. The organism can easily enter the food processing chain and cause contamination of food products or equipment. *B. cereus* produces emetic and diarrheagenic enterotoxins, which can cause vomiting and diarrhea when consumed. Consequently, routine testing for *B. cereus* in food processing or handling areas is recommended. Compact Dry™ BC can be used for the quantitative detection of *Bacillus cereus* in foods or other samples in place of MYP Agar.⁽¹⁻⁵⁾

Compact Dry™ BC is a ready-to-use chromogenic medium for performing *Bacillus cereus* colony counts. The product contains dehydrated culture media and a cold water-soluble gelling agent in a non-woven cloth matrix. The medium is instantly hydrated when inoculated with a sample, and capillary action diffuses the sample evenly over the matrix to form a gel within seconds. Compact Dry™ BC contains a chromogenic substrate that yields a colored reaction when utilized by cells and permits the differentiation of *B. cereus* from background flora in the sample.

Compact Dry™ BC has been validated to ISO 16140 against the reference method (EN ISO 7932:2004) for the detection of *Bacillus cereus* in foods (e.g. meat products, fish and seafood, milk products, pasta, and bakery items). Compact Dry™ BC is also MicroVal (no. 2019LR87) and NordVal (no. 045) validated and can be routinely applied by certified users. The ready-to-use trays save space and greatly reduce the time needed to perform microbiological testing. Compared to other commonly used culture systems, Compact Dry™ has a longer shelf life, can be stored at room temperature, does not require manual sample spreading, is rigid, stackable and easy to label, and allows for direct colony picking for further subculture.

FORMULA

Compact Dry™ BC contains dehydrated culture media, a gelling agent, and a chromogenic substrate to assist in differentiating *Bacillus cereus* from foods or other samples.

Final pH 7.0 +/- 0.2 at 25°C.

STORAGE AND SHELF LIFE

Storage: Upon receipt, store at 1-30°C. away from direct light. Media should not be used if there are any signs of deterioration, contamination, or if the expiration date has passed. Product is light and temperature sensitive; protect from light, excessive heat, moisture, and freezing. If foil pouch is opened and not all plates are used, return remaining plates to pouch and reseal until next use. Opened packages should be used as soon as possible.

The expiration dating on the product label applies to the product in its intact packaging when stored as directed. The product may be used and tested up to the expiration date on the product label and incubated for the recommended quality control incubation times.

Refer to the document "[Storage](#)" for more information.

PRECAUTIONS

This product may contain components of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not guarantee the absence of transmissible pathogenic agents. Therefore, it is recommended that these products be treated as potentially infectious, and handle observing the usual universal blood precautions. Do not ingest, inhale, or allow to come into contact with skin.

This product is for laboratory use only. It is to be used only by adequately trained and qualified laboratory personnel. Observe approved biohazard precautions and aseptic techniques. All laboratory specimens should be considered infectious and handled according to "standard precautions." The "Guidelines for Isolation Precautions" is available from the Centers for Disease Control and Prevention at www.cdc.gov/ncidod/dhqp/gl_isolation.html.

For additional information regarding specific precautions for the prevention of the transmission of all infectious agents from laboratory instruments and materials, and for recommendations for the management of exposure to infectious disease, refer to CLSI document M-29: *Protection of Laboratory Workers from Occupationally Acquired Infections: Approved Guideline*.

Sterilize all biohazard waste before disposal.

Refer to the document "[Precautions When Using Media](#)" for more information.

Refer to the document [SDS Search](#) instructions on the Hardy Diagnostics' website for more information.

PROCEDURE

Prior to Use: Refer to listed references for appropriate methods of collection, preparation, and dilution of specimens under investigation. ⁽¹⁻⁵⁾

For environmental samples or to test uneven surfaces of equipment, swab surface of test area with [EnviroTrans™](#) (e.g. Cat. no. SRK05 or SRK35).

For raw material and food testing, prepare and dilute samples using an appropriate diluent such as [Dilu-Lok™ II](#) (e.g. Cat. no. D590 or D599).

General Dilution Guidelines:

For Making 1:10 Serial Dilutions

1. Using a sterile pipet or scoop, aliquot 10ml or 10gm of test sample to a 90ml pre-filled Dilu-Lok II™ dilution vial to yield a 1:10 dilution. Mix thoroughly.
2. From the 1:10 dilution vial in step 1, use a fresh sterile pipet and aliquot 10ml from this dilution vial into a second 90ml pre-filled Dilu-Lok II™ vial to yield a 1:100 dilution. Mix thoroughly.
3. Continue aliquoting 10ml dilutions into 90ml pre-filled Dilu-Lok II™ vials until the desired concentration of test sample is achieved. Each subsequent dilution increases by a factor of 10. A separate sterile pipet should be used with each dilution.

For Making 1:100 Serial Dilutions

1. Using a sterile pipet or scoop, aliquot 1ml or 1gm of test sample to a 99ml pre-filled Dilu-Lok II™ dilution vial to yield a 1:100 dilution. Mix thoroughly.
2. From the 1:100 dilution vial in step 1, use a fresh sterile pipet and aliquot 1ml from this dilution vial into a second 99ml pre-filled Dilu-Lok II™ vial to yield a 1:10,000 dilution. Mix thoroughly.
3. Continue aliquoting 1ml dilutions into 99ml pre-filled Dilu-Lok II™ vials until the desired concentration is achieved. Each subsequent dilution increases by a factor of 100. A separate sterile pipet should be used with each dilution.

Method of Use Direct Inoculation:

1. Remove the set of four trays from the foil pouch and separate each individual tray by gently bending along the connecting edge until each tray snaps free. Alternatively, if setting up a dilution series of the same sample, trays can be left connected to facilitate reading similar samples. Trays that are not used immediately should be resealed in the foil pouch. Refer to the "Storage and Shelf Life" section for proper storage of unused trays.
2. Remove the lid of the tray using two fingers to hold down one end of the lid and the thumb to lift the opposite end. Lids are easier to remove using a "peel back" method as opposed to a "pull off" method.
3. Inoculate by pipetting 1ml of sample directly to the center of a dry tray well, being careful not to touch the surface of the matrix with the pipet tip. Once dispensed, the sample will automatically diffuse across the surface by capillary action to form a gel; manual spreading of the inoculum is discouraged. Remember to account for the sample inoculum when calculating the dilution series.
4. Replace the lid and label the tray with appropriate information, including the sample dilution factor.
5. Invert the tray and incubate, upside down with the medium on top, at 30° +/- 1°C for 24 +/- 2 hours. NOTE: Use the appropriate temperature/time designation according to the legal specification of the prescribed food analysis regulation.
6. Count colonies using the Hardy Diagnostics Wizard™ Compact Dry™ plate reader ([Cat. no. CDR1](#)) designed exclusively for use with Compact Dry™. See the Wizard™ Compact Dry™ instruction manual. Alternatively, colonies can be counted when illuminated from the backside of the tray to calculate CFU/ml using the Scan® 100 colony counter ([Cat. no. 435000](#)) or comparable back lighting. If the colony count is high, use the 1cm x 1cm molded grid on the back of the tray to assist in colony counting. Use a sheet of white paper with gridded lines to diffuse the light if the molded grids in the tray are difficult to visualize with a light box.

Method of Use Membrane Filtration:

1. Remove the lid on a Compact Dry™ plate and pipette 1.0ml of sterile purified water into the middle of the plate to activate the matrix just before use.
2. Using a membrane filtration set-up, filter 100-250ml of a water sample under reduced pressure through a sterile 0.45µm membrane filter.
3. Keep the set-up running and rinse the inside of the filter funnel using 20-30ml of sterile purified water. Repeat this step two to three times to ensure all of the sample has been filtered through the membrane.
4. Remove the membrane from the funnel using sterile forceps (e.g. [Cat. no. 800000](#)) and apply it filter side up (trap side away from the matrix) to the surface of the pre-moistened Compact Dry™ plate. Gently press the membrane onto the surface, making sure to remove bubbles so the membrane is completely flush and centered on the matrix.
5. Replace the lid and incubate the tray right-side-up using the information outlined above.
6. Count colonies and record results as outlined above.

INTERPRETATION OF RESULTS

After incubation, read trays using the Wizard™ Compact Dry™ plate reader ([Cat. no. CDR1](#)) or read colonies against a white or illuminated background such as with the Scan® 100 colony counter (Cat. no. 435000) or comparable back lighting.

Bacillus cereus forms light blue to blue colonies on the medium. Bacteria other than *B. cereus* may grow and produce white colonies. Only light blue to blue colonies should be counted. *Bacillus thuringiensis* or closely related species with similar biochemical properties may grow and form light blue to blue colonies. See the Limitations section for more information. For methods requiring a confirmation step, subculture light blue to blue colonies to a Blood Agar Plate with 5% Sheep Blood ([Cat. no. A10](#)) for further testing.

LIMITATIONS

During inoculation, do not touch the surface of medium and be careful to avoid any contamination by airborne microorganisms.

During incubation, keep cap tight on plates to avoid any possible dehydration.

A dilution may be needed when the sample has a dark color.

When the sample is viscous (thick), pipetting the sample on several points on a plate or an additional dilution may be needed for an even suspension.

When the sample contains an enzyme, it may react with the enzyme substrate in the dry sheet and affect the color.

If the nature of sample does affect the reaction of the medium, inoculate only after the factor is eliminated by means of dilution and other techniques. (e.g. samples with high viscosity, colored, reactive with chromogenic substrate, and with a high or low pH).

It is recommended to use a stomacher and filter homogenized sample afterwards to eliminate carry-over of tiny particles of foodstuff onto the surface of the medium.

Counting colonies may be difficult against a dark background. For best results, count colonies using the Hardy Diagnostics Wizard™ Compact Dry™ plate reader ([Cat. no. CDR1](#)) or with the tray held against a white or illuminated background such as with the Scan® 100 colony counter (Cat. no. 435000).

If using a light box, the molded grid lines or colonies may be difficult to view due to excessive brightness. Diffuse the light using a sheet of white, gridded (1cm x 1cm) paper underneath the tray to facilitate colony counting.

Colonies are not distinguishable on trays if concentrations are above 100 CFU/ml, as high colony counts will result in the whole surface becoming colored. The sample should be diluted to a concentration of less than 100 CFU/ml for best use.

Failure to pre-moisten the Compact Dry™ plate matrix prior to performing membrane filtration will result in poor performance.

Placing the filter (i.e. trap) side down on the matrix when performing membrane filtration may result in poor performance and an inability to count individual colony forming units (CFUs).

It is recommended the membrane filtration method employed by the user be validated prior to use with Compact Dry™.

If present in the sample, closely related species such as *B. thuringiensis*, *B. weihenstephanensis*, and *B. mycoides* may also grow and produce light blue to blue colored colonies.

Refer to the document "[Limitations of Procedures and Warranty](#)" for more information.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as loops, swabs such as EnviroTrans™, applicator sticks, scoops, dilution buffers such as Dilu-Lok™ II, other culture media, Wizard™ Compact Dry™ plate reader ([Cat. no. CDR1](#)), Scan® 100 colony counter (Cat. no. 435000), incinerators, forceps ([Cat. no. 800000](#)), membrane filtration laboratory products and equipment, and incubators, etc., as well as serological and biochemical reagents, are not provided.

QUALITY CONTROL

End users should anticipate the following typical performance characteristics when testing with Compact Dry™.

Test Organisms	Inoculation Method*	Incubation			Results
		Time	Temperature	Atmosphere	
<i>Bacillus cereus</i> ATCC® 13061	J	24hr	30°C	Aerobic	Growth; light blue to blue colonies
<i>Bacillus thuringiensis</i> ATCC® 33679	J	24hr	30°C	Aerobic	Growth; light blue to blue colonies
<i>Staphylococcus aureus</i> ATCC® 29213	J	24hr	30°C	Aerobic	Growth; white colonies
<i>Bacillus subtilis</i> ATCC® 6633	J	24hr	30°C	Aerobic	Inhibited

* Refer to the document "[Inoculation Procedures for Media QC](#)" for more information.

USER QUALITY CONTROL

Check for signs of contamination and deterioration. Users of commercially prepared culture media may be required to perform quality control testing to demonstrate growth or a positive reaction and to demonstrate inhibition or a negative reaction (where applicable). See the following reference for more specific information. ⁽¹⁻⁵⁾

PHYSICAL APPEARANCE

Compact Dry™ BC should appear dry, free of particles, and off white in color.

REFERENCES

1. Association of Official Analytical Chemists. 2005. *Official Methods of Analysis*, 18th ed. AOAC, Washington, D.C.
2. American Public Health Association. 2012. *Standard Methods for the Examination of Dairy Products*, 22nd ed. APHA, Washington, D.C.
3. APHA Technical Committee on Microbiological Methods for Foods. 2001. *Compendium of Methods for the Microbiological Examination of Foods*, 4th ed. APHA, Washington, D.C.
4. U.S. Food and Drug Administration. *Bacteriological Analytical Manual*. AOAC, Arlington, VA.
www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm2006949.htm
5. Teramura, H., M. Uchida, and H. Kodaka. 2011. CompactDry X-BC for the Enumeration of *Bacillus cereus* in Food Samples. *Biocont Sci.* 16:2, 73-77.

ATCC is a registered trademark of the American Type Culture Collection.
Compact Dry is a trademark of Nissui Pharmaceutical Co., Ltd.
Scan is a registered trademark of Interscience for Microbiology.

MicroVal approval no. 2019LR87
NordVal approval no. 045

IFU-000811[A]

Manufactured for and distributed by:



1430 West McCoy Lane, Santa Maria, CA 93455, USA

Phone: (805) 346-2766 ext. 5658

Fax: (805) 346-2760

Website: www.HardyDiagnostics.com

Email: TechService@HardyDiagnostics.com

[Ordering Information](#)

Distribution Centers:

California · Washington · Utah · Arizona · Texas · Ohio · New York · Florida · North Carolina

The Hardy Diagnostics manufacturing facility and quality management system is certified to ISO 13485.

Copyright© 2020 by Hardy Diagnostics. All rights reserved.

HDQA 22071 [E]