InSite Listeria

Environmental Listeria Species Test

Part No: IL100 (100 tests), IL050 (50 tests)

Description/ Intended Use:

InSite Listeria is a screening test for *Listeria* spp., intended for use on food contact surfaces and food processing equipment after cleaning to detect the presence of *Listeria* species. A color change of the media from yellow/amber to light brown/black is considered presumptive positive.

Principle:

InSite Listeria contains a patented formula of antibiotics, growth enhancers and color changing compounds. The antibiotics inhibit most non-*Listeria* microorganisms while growth enhancers provide recovery nutrients to support growth of sub-lethally injured *Listeria*. Indicator compounds turn broth from yellow to black by utilizing β -glucosidase enzyme produced by *Listeria* species. A brown/black color after 24 – 48 hours at 37 °C indicates a presumptive positive result for *Listeria* spp.

Required Materials (Not Provided):

- Incubator set to 37 ± 1 °C







Directions:

1. When collecting sample, make sure to use aseptic technique. Do not touch swab or inside of sampling device. Holding swab tube firmly, twist and pull top of swab out of swab tube. Foam tip swab is pre-moistened; condensation may be visible on inside of swab tube - this is normal. Thoroughly swab a standard 30 x 30 cm (12 x 12 inch) area of interest for a typical flat surface. Rotate swab as sample is being collected to ensure maximum sample pickup and apply sufficient pressure to create flex in swab shaft. For irregular surfaces, ensure swabbing technique remains consistent for each test and swab a large enough area to collect a representative sample.
2. After swabbing desired test area, place swab back in swab tube and close device firmly.
3. To activate InSite, hold swab tube firmly and use thumb and forefinger to break Snap-Valve by bending bulb forward and backward. Squeeze bulb 3 to 4 times, expelling all liquid down into tube.
4. Gently massage bottom of tube by squeezing tube 3 times, then shake for 3 seconds. This will help release cells from swab and displace air bubbles.
5. Incubate for 24 – 48 hours at 37 ± 1 °C. Observe color change and record result. Presumptive positive result can be read as early as 24 hours. Results cannot be considered negative until sample has been incubated for 48 hours.

Interpretation of Results:

- If media color changes from yellow/amber to brown/black within 24 – 48 hours of incubation at 37 ± 1 °C, sample is presumptive positive for *Listeria* spp.
- A slight darkening of broth could suggest broth is beginning to change and may require an additional 4 hours of incubation. As with all laboratory protocols, use of a negative control is recommended.
- Yellow / amber media with black spots on swab bud should be considered presumptive negative for *Listeria* spp.

Note: A naturally occurring orange precipitation may be seen in some devices. This is due to excess iron falling out of solution. To alleviate this, device can be shaken gently to dissolve, either before or after incubation, and used as normal. This will not affect test device performance.

Confirmation:

Presumptive positive samples can be confirmed by streaking sample onto commonly used selective *Listeria* agar plates such as Modified Oxford agar, Palcam agar, or any other recognized confirmatory procedure. Typical *Listeria* colonies on selective agar plates could then be further analyzed by more definitive tests such as microscopy, biochemical tests, etc.

Storage & Shelf Life:

- Store devices at 2 – 8 °C (35 – 46 °F)
- Devices have a 12 month shelf life.
- Check expiration date on label.

Disposal:

Disinfect before disposal. InSite devices can be disinfected by autoclaving, incinerating, or by soaking in 20% bleach for 1 hour. Then, they can be placed in the trash. Alternatively, InSite devices may be discarded at a biohazard waste disposal facility.

Safety & Precautions:

Components of InSite devices do not pose any health risk when used correctly. Used devices confirming positive results may be a biohazard and should be disposed of safely in compliance with Good Laboratory Practice and Health and Safety Regulations.

- InSite is intended to be used on product contact surfaces after cleaning. Heavily contaminated surfaces such as drains and floors will give false positive results.
- *Listeria monocytogenes* (*L. mono*) is a dangerous human pathogen. When handling samples that possibly contain *L. mono*, extreme care should be taken to contain samples and enriched samples (presumptive positive InSite devices). Immuno-compromised individuals and pregnant women are particularly endangered by exposure to *L. mono* and should not be allowed in the vicinity of testing.

Hygiena Liability:

As with any culture medium, InSite Listeria results do not constitute a guarantee of quality of food, beverage products or processes that are tested with these devices. Hygiena will not be liable to user or others for any loss or damage, whether direct or indirect, incidental or consequential from use of these devices. If this product is proven to be defective, Hygiena's sole obligation will be to replace product, or at its discretion, refund the purchase price. Promptly notify Hygiena within 5 days of discovery of any suspected defect and return product to Hygiena. Please call Customer Service for a Returned Goods Authorization Number.

Contact Information:

If more information is required, please visit us at www.hygiena.com or contact us at:

Hygiena - Americas
Phone: 805.388.8007
Fax: 805.388.5531
Email: info@hygiena.com

Hygiena - International
Phone: +44 (0)1923 818821
Fax: +44 (0)1923 818825
Email: enquiries@hygiena.com

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