



HARDY DIAGNOSTICS

# Dilu-Lok II™

The Dilu-Lok II™ eliminates the time and money spent on preparing your own dilution vials! Our pre-filled dilution vial is designed especially for the food, dairy, pharmaceutical, and water industries. The vials are filled to 25ml, 90ml, or 99ml with a variety of buffering and enrichment solutions. Each lot is meticulously tested for pH, sterility, fill volume, and toxicity.

**NEW!** Call to receive a FREE sample pack to evaluate in your lab!

## Save time and money with the Dilu-Lok II™ pre-filled dilution vial



### Simple...

Easy-to-remove, perforated, shrink wrap ensures product integrity.



### One-Handed...

Hinged lid allows for easy opening.



### Large Opening...

A wide 45mm opening accommodates most sample sizes.



### Leak-Proof...

Recently re-engineered for an improved seal.

Please order in multiples of 50.

Product	Cat. no.
Butterfield's Phosphate, 99ml, 50 per box .....	D599
Butterfield's Phosphate, 90ml, 50 per box .....	D590
Butterfield's Phosphate, 25ml, 50 per box .....	D125
Peptone Water 0.1%, 99ml, 50 per box .....	D299
Peptone Water 0.1%, 90ml, 50 per box .....	D290
Phosphate Buffer with MgCl, 99ml, 50 per box ...	D699
Phosphate Buffer with MgCl, 90ml, 50 per box ...	D690



Santa Maria, CA | Salt Lake City, UT | Phoenix, AZ | Tacoma, WA  
Lima, OH | Dallas, TX  
Phone: 800.266.2222  
Fax: 805.346.2760  
Website: [www.HardyDiagnostics.com](http://www.HardyDiagnostics.com)  
Email: [Sales@HardyDiagnostics.com](mailto:Sales@HardyDiagnostics.com)