

Sure Tight™ Sample Cup

*Tight Sealing Lid
Ensures that Your Sample will
Arrive in the Lab Safe and Dry!*



The SureTight™ from Hardy are tough, leak proof, sterile, sampling cups.

Your Sample Cup "Dream Come True"

SureTight™ represents a new advance in container technology.

The patented lid is designed to securely lock in the sample. When you use the SureTight™, you can be confident that your sample is safely contained in the cup and will remain that way throughout transport to the lab!

- ◆ Tamper-proof sterility tab ensures sample integrity.
- ◆ The cups feature a wide 53mm opening for larger samples.
- ◆ SureTights™ are constructed of tough polypropylene with thick walls making them virtually indestructible.
- ◆ Ergonomic design features lids with ridges that allow the technician to easily open the cups.
- ◆ A paper label is attached for easy sample identification. Graduation marks on the side allow for volume measurement.
- ◆ SureTights™ are stackable and interlocking for efficient use of space.
- ◆ Economical - By far, SureTight™ represents the best value!

90ml (3oz)	Sterile	400/case	Cat# 243810
------------	---------	----------	-------------



Tel: (800) 266-2222 Website: HardyDiagnostics.com
Fax: (805) 346-2760 Email: Sales@HardyDiagnostics.com

Santa Maria, CA Tacoma, WA Salt lake City, UT Phoenix, AZ Lima, OH