

# Evaluation of StrepB Carrot Broth versus Lim Broth for Detection of Group B Streptococcus Colonization Status of Near-Term Pregnant Women<sup>▽</sup>

D. L. Church,<sup>1,2,3\*</sup> Heather Baxter,<sup>4</sup> Tracie Lloyd,<sup>1</sup> Beverley Miller,<sup>1</sup> and Sameer Elsayed<sup>3</sup>

*Calgary Laboratory Services<sup>1</sup> and Departments of Pathology & Medicine,<sup>2</sup> Medicine,<sup>3</sup> and Family Medicine,<sup>4</sup> Calgary Health Region and the University of Calgary, Calgary, Alta., Canada*

Received 24 March 2008/Returned for modification 6 May 2008/Accepted 27 May 2008

**The performance of StrepB Carrot Broth (SCB) versus group B Lim broth (LIM) for detection of group B streptococcus (GBS) colonization status in near-term pregnant women (35 to 37 weeks of gestation) was evaluated. Dually collected vaginal/rectal swabs from 279 women enrolled from a single large maternity clinic were analyzed. Fifty (18%) women were colonized by GBS according to both methods. SCB had excellent diagnostic performance compared to LIM, with sensitivity, specificity, positive predictive value, and negative predictive value of 92%, 100%, 100%, and 98.3%, respectively. Improved diagnostic efficiency due to direct reporting of GBS cases based on an orange color change in the SCB decreased overall labor and material costs.**