



Sterility Testing Media

<USP 71>

Media is tested to USP <71> standards.
 Media is manufactured in an ISO 13485 certified facility.
 Ideal for use with the Sterisart® Sterility Testing System.

Tryptic Soy Broth, USP

Recommended for use as a general purpose medium for the isolation and cultivation of a wide variety of bacteria and fungi.
 100ml glass bottle with needle port septum, 100ml fill, 20/pk.....U46



Fluid Thioglycollate Broth Medium with Indicator, USP

For the cultivation of microorganisms and USP sterility procedures.
 100ml glass bottle with needle port septum, 100ml fill, 20/pk.....U84



Go to HardyDiagnostics.com for a complete listing of bottle and jar media. Custom containers and fill volumes are available upon request.

Fluid A, USP

A rinsing and diluting fluid for use in USP testing protocols. Contains peptone 0.1%.
 100ml glass bottle with needle port septum, 100ml fill, 20/pk.....U109



Fluid D, USP

A rinsing and diluting fluid for use in USP testing protocols. pH 7.1. Contains peptic digest of animal tissue and polysorbate 80.
 100ml glass bottle with needle port septum, 100ml fill, 20/pk.....U110



1430 West McCoy Lane,
 Santa Maria, CA 93455
 800-266-2222
Sales@HardyDiagnostics.com