

QUICKSLIDE™

Automated Gram Stainers

A Division of
HARDY
DIAGNOSTICS



 Online Ordering Available



QUICKSLIDE™ GRAMPRO 1

Cat. no AGS-1000

ELECTRO-OPTICAL PROCESS

The GramPRO 1™ automatically adjusts for sample thickness, has individual preference controls for color intensity and is operator friendly.

FAST, HANDS-FREE STAINING

Totally hands free, each slide is automatically stained, decolorized and counter stained in 3.5 to 4.5 minutes. The GramPRO 1™ has an automatic cleaning cycle, has internal diagnostic software.

SMALL COMPACT SIZE

The GramPRO 1™ is small enough to fit on any lab counter top.

MINIMAL DAILY MAINTENANCE

Easy and fast. It's as simple as the push of a button.

WHAT ARE THE BENEFITS OF THE GramPRO 1™ UNIT OVER MANUAL GRAM STAINING?

- The GramPRO 1™ is not "technique dependent". You get consistency no matter who processes the slide. With our patented electronic eye, decolorization is perfect every time.



FEATURES

- Patented Decolorization process that optically monitors run-off fluid at a rate of 20 times per second
 - Totally hands free process – Soft touch keypad
 - Stain times from 3.5 minutes to 4.5 minutes
 - Adjustable stain settings for color intensity
 - Automatic calibration of each sample
 - Adjustable sleep mode for energy efficiency
 - Stain resistant metal cabinet fits on counter top
 - Consistency for all sample types
 - Automatic cleaning cycle
 - Internal diagnostic software
 - Audible on/off alarm to signal end of test
 - Decolorize only cycle
 - Manufactured in USA
-
- There is no cross-contamination. The GramPRO 1™ uses fresh stain for each sample slide.
 - No more filtering of stains as with the manual process. This saves time and eliminates the mess.
 - The GramPRO 1™ has independently adjustable stain times on both the long and short cycles. It also has the option of decolorizing only if the operator chooses to do so.
 - The GramPRO 1™ allows the user to do other tasks while the slide is processing.
 - No technical expertise is required to do a quality Gram Stain. As long as the sample is properly heat fixed, anyone can perform the Gram Staining process with confidence.
 - Fewer fluids are used saving the lab disposal fees. The GramPRO 1™ does not have to be placed by sink – it can be placed beside a microscope.
 - Spend more time testing and less time cleaning your machines. Check out the competitors' recommended scheduled maintenance and you will appreciate the savings being offered by the GramPRO 1™.

TOP SELLER!





Cat. no QS-001

The HemaPRO™ Slide Stainer represents a significant advancement in automated slide staining and is a valuable addition to any laboratory, whether it is the primary stainer used by small hospitals, physicians labs, stat labs, or used as a stat/backup instrument in larger labs. The HemaPRO™ can also be utilized in veterinarian labs for a variety of biological specimens.

NEW FEATURE!

In addition to the soft agitation before the Rinse is drained, HemaPRO™ has an added soft agitation feature that occurs before the Buffer is drained. Through extensive testing we have determined this added feature puts off a cleaner slide with the same great staining characteristics of our modified Wright Giemsa stain formula. The unit needs no daily maintenance and only requires a tubing change once or twice annually.



FEATURES

- State-of-the-art microprocessor, increased speed and ease of operator use
- Operator prompt LCD display
- Soft touch color-coded keypad
- Adjustable stain and buffer timing – so important for consistent lab results with multiple users
- Adjustable fill level
- “Cycles Remaining” display – computer monitors reagent and prompts when to reorder
- Injection molded case with built-in reagent platform
- Stainless steel reagent cannulas – oxidation resistant, increase system durability
- Color-coded reagent tubing for ease of operation
- Special formulated reagent kit – 125 tests per kit in one self-contained package

SMALL COMPACT SIZE

The HemaPRO™ is one of the smallest automated slide stainers available. It requires only one square foot of counter space. In addition to its small size, the HemaPRO™ reagent system is included in the one square foot.

FAST STAINING

The HemaPRO™ can stain a blood smear in as little as 60 seconds. In addition to the rapid staining, the HemaPRO™ eliminates cross contamination seen in the dip method and yields consistent results no matter who stains the slide.

PERFORMS MULTIPLE STAINING PROCEDURES

In addition to staining peripheral blood smears, the HemaPRO™ has a special program that allows it to stain bone marrow samples. The HemaPRO also allows operator to adjust both stain and buffer times.



1430 West McCoy Lane
Santa Maria, CA 93455

