

Our pre-filled dilution vials are designed especially for the food, dairy, pharmaceutical, and water industries. The 145ml vials have a hinged cap and are filled to 25ml, 90ml, or 99ml with a variety of buffering and enrichment solutions. Each lot is meticulously tested for pH, sterility, and fill volume.

*Save time
and money
spent
preparing
your own
dilutions!*

Dilu-Lok™ Dilution Vials




**DILU-
LOK™**

Prefilled Dilution Vials

Butterfield's Phosphate Buffer, 25ml, 50/pk	D125
Butterfield's Phosphate Buffer, 90ml fill, 50/pk	D590
Butterfield's Phosphate Buffer, 99ml fill, 50/pk	D599
Peptone Water 0.1%, 90ml, irradiated, 50/pk	D290
Peptone Water 0.1%, 99ml, irradiated, 50/pk	D299
Phosphate Buffer with Magnesium Chloride, 90ml fill, 50/pk	D690
Phosphate Buffer with Magnesium Chloride, 99ml fill, 50/pk	D699
TAT (Tryptone Azolectin Tween®) Broth, 90ml, irradiated, 50/pk	D077
Buffered Peptone Water, 90ml fill, irradiated, 50/pk	D080
Buffered Peptone Water, 99ml fill, irradiated, 50/pk	D089
Buffered Peptone Water with 1% Tween® 20, 90ml fill, irradiated, 50/pk	D085
Deionized Water, 50ml fill, 50/pk	D050
Deionized Sterile Water, 90ml fill, 50/pk	D090
Deionized Sterile Water, 99ml fill, 50/pk	D099

- **Simple** - Easy-to-remove, perforated, shrink seal to ensure product integrity
- **Large Opening** - A wide 45mm opening to accommodate most sample sizes
- **Certified** - Each case delivered with a certificate of quality
- **One-Handed** - Hinged lid for easy opening
- **Leak-Proof** - New improved seal



Easy Opening Flip-top

Convenient Write-on area

Wide mouth opening

Hinged lid

99ml

90ml

Clarified Polypropylene to easily visualize the sample

Recyclable plastic

HARDY
DIAGNOSTICS
A Culture of Service™

100%
Employee Owned

ESTABLISHED
1980

Made in the USA

ISO 13485
Medical Devices
Quality Management
FM 572526

Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.

Headquarters
1430 West McCoy Lane
Santa Maria, CA 93455
800.266.2222
Sales@HardyDiagnostics.com
HardyDiagnostics.com

Distribution Centers
Santa Maria, California
Olympia, Washington
Salt Lake City, Utah
Phoenix, Arizona
Dallas, Texas
Springboro, Ohio
Lake City, Florida
Albany, New York
Raleigh, North Carolina

Copyright © 2019 Hardy Diagnostics