


HARDY
DIAGNOSTICS
A Culture of Service™



TC (Total Count)
1) Total Count = 49



Total Count = 49
Unknown Count = 0



The Wizard
COMPACT DRY READER

The Wizard
CompactDry™ Reader



Automate your colony counts and expand your throughput capabilities.

Your choice of data retrieval:

- Printout
- Thumb Drive
- LIMS Connection



Cat. no. CDR1

Improve your work flow!

- Create **consistency** from tech to tech with automation your team can rely on.
- **Increase** daily **throughput** by automating the tedious task of colony counting.
- **Standardize** your testing procedures.

For use with these CompactDry™ plates:



E. coli and Coliforms
Cat. no. EC100



Total Count
Cat. no. TPC100



Yeast and Mold
Cat. no. YM100



Yeast and Mold Rapid
Cat. no. YMR100



Enterobacteriaceae
Cat. no. ETB100



Listeria
Cat. no. LS100



Staphylococcus aureus
Cat. no. XSA100

Simplify data storage, sample entry, and reporting with the accessories below.

Accessories



USB Wireless Keyboard

Includes
USB dongle.
Cat. no. 120103



Barcode Scanner

Includes USB Connector
and stand.
Cat. no. 120102



Thermal Printer

Includes USB Connector
and power supply.
Cat. no. 120101



Contact us:

800.266.2222

Sales@HardyDiagnostics.com

HardyDiagnostics.com

1430 West McCoy Lane
Santa Maria, CA 93455

Distribution Centers

Santa Maria, California
Olympia, Washington
Salt Lake City, Utah
Phoenix, Arizona
Dallas, Texas
Springboro, Ohio
Raleigh, North Carolina
Lake City, Florida
Albany, New York



100%
Employee Owned



Made in the USA

