



MANUFACTORY
PURVEYORS OF INDUSTRY
CALIFORNIA

2421 S. WESTGATE ROAD
SANTA MARIA, CA 93455
T: 805.354.7511
F: 805.922.4980

PRODUCT DATA SHEET

MD3033.1 Face Shield

I. OVERVIEW



Well suited for protection against splashes

Provides clear, distortion-free view with anti-fog coating on inside and outside of shield

Provides complete coverage including sides of face, eyes, nose and mouth

Includes a comfortable, skin friendly forehead cushion

Fits easily over goggles or glasses and masks

One size fits all due to adjustable head strap

II. INFORMATION

Product Name:	Disposable Face Shield
Item no.:	MD3033.1
Product Type:	General Use
Units per pack:	10
Unit weight:	28-29grams
Product Dimensions:	Shield: 330 x 220 x 0.20 mm Foam: 260 x 30 x 30 mm Elastic: 340 x 150 mm
Imported:	Yes

III. SPECIFICATIONS

Materials:	Shield: Food-grade PET with removable protective film Elastic: Polyester rubber Foam: High dense, open foam
Characteristics:	Anti-fog treatment Fiberglass-free Latex-free Lightweight Distortion-free Foam headband Optically clear

IV. CONTACT INFORMATION

Company Name:	The Manufactory, USA
Address:	2421 Westgate road, Santa Maria, California 93455, USA
Contact no:	(805) 922.2250
Website:	www.manufactoryusa.com

Cleaning instructions: When disinfecting or sanitizing our face shields, the following products should be used: Steriplex SD, Safetec Surface Wipes, Perodox RTU, Prostat Sterile Wipes, hydrogen peroxide (pharmacy grade), bleach (household grade), or isopropyl alcohol (diluted with 30% water).

Issued Date: 04.15.2020