

anaero₂GRO™

Pre-Reduced Anaerobic Culture Media

Features	AnaeroGRO™	Anaerobe Systems
Performance	Superior	
Anaerobically packaged	Yes	Yes
Oxygen scavenger included	Yes	No (Peice of paper only)
Desiccant included	Yes	No
All formulas offered	Yes	Yes
Primary set-up packs (multi-dish) offered	Yes	Yes
Pre-reduced	Yes	Yes
Anaerobically poured	No	Yes
Easy tear-notch on pouch	Yes	No
Manufactured in an ISO certified facility	Yes	No
Room temperature storage	Yes	Yes
Shelf life 90 days for plates	Yes	Yes
Can be ordered with standard media	Yes	No
Number of distribution centers in the US	9	1

*Culture media that is exposed to environmental oxygen leads to a build-up of reactive oxygen species (ROS) that initiate damaging free radical reactions, which inhibit the growth of anaerobic bacteria. Therefore, ingredients have been added to the AnaeroGRO media to neutralize the growth inhibiting toxic effects of peroxide and other reactive oxygen species (ROS) that may develop during the medium's brief exposure to oxygen after it is sterilized and before it is packaged in an oxygen free environment.



060619sg



FM 572526

Hardy Diagnostics has a Quality Management System that is certified to ISO 13485 and is a FDA licensed medical device manufacturer.



Made in the USA



100%
Employee Owned

Headquarters

1430 West McCoy Lane
Santa Maria, CA 93455
800.266.2222
Sales@HardyDiagnostics.com
HardyDiagnostics.com

Distribution Centers

Santa Maria, California
Olympia, Washington
Salt Lake City, Utah
Phoenix, Arizona
Dallas, Texas
Springboro, Ohio
Lake City, Florida
Albany, New York
Raleigh, North Carolina

Copyright © 2019 Hardy Diagnostics