

# QUICKSLIDE™

## Quick Start Guide

### GramPRO 1 Automated Gram Stainer

**Congratulations!** You have made an excellent choice for your lab. Thank you for partnering with **Hardy Diagnostics**. Now, let's get started so you can begin using your new **GramPRO 1**. Just follow these simple instructions....

### You should have received

GramPRO 1 \* Quick Start Thumb Drive \* User Manual \* Power Cord



### First things first

- 1. Plug the Quick Start thumb drive into your computer's USB port** - Open thumb drive titled "GramPRO 1."
- 2. Register your GramPRO 1 warranty** - We provide a free one-year standard warranty which begins on the date of purchase. To activate your standard warranty open the "Register GramPRO 1" link in the Quick Start thumb drive or visit [HardyDiagnostics.com/Register-GramPRO](http://HardyDiagnostics.com/Register-GramPRO).
- 3. Get to know your GramPRO 1** - Open the "GramPRO 1 Videos" link to review the instructional videos on YouTube.
  - **Installation** - Review the "GramPRO 1 Installation" video to start up your GramPRO 1 or reference the "First Time Instrument Set-Up" section in your User Manual.
  - **Operation** - Review the "GramPRO 1 Operation" video to learn how to operate your GramPRO 1. Operation instructions are also detailed in the "First Time Instrument Set-Up" section of your User Manual.
  - **Maintenance** - Review the "GramPRO 1 Maintenance" video to learn how to maintain your GramPRO 1. Maintenance instructions are also available in your User Manual.
  - **GramPRO 1 Website** - If you would like more information about the GramPRO 1 at any time, please open the GramPRO 1 website link or visit [HardyDiagnostics.com/GramPRO](http://HardyDiagnostics.com/GramPRO).

**If you would like phone support while setting up your instrument, please call our Technical Support team at 800-266-2222 ext. 5598.**

This Quick Start Guide is not intended to replace the User Manual. Please visit our website at [HardyDiagnostics.com](http://HardyDiagnostics.com) to obtain a copy of the most recent User Manual.

# QUICKSLIDE™

## Supplementary Products



**GP1SK**  
GramPRO 1 Stain kit:  
Required to commence GP1 stain cycle



**PB72**  
PROBond slides:  
For superior specimen adhesion



**GP1RT**  
GramPRO 1 tubing kit:  
Change tubing every six months or 1,100 cycles



**QS1016**  
Methanol:  
For specimen fixation



**QS2000**  
Gram's Iodine:  
To disinfect rinse water



**QS1001**  
Foam tipped swabs:  
Used for cleaning the cuvette of stain residue



**QS4000**  
System cleaner:  
This is required for the GramPRO 1 scrub cycle which must be performed every 15 days or 300 cycles



**QS0700**  
QSlide PROBond QC slides:  
Control slides with Gram negative and positive bacteria

For more information on the GramPRO 1 and supplementary products, please visit [HardyDiagnostics.com/GramPRO](http://HardyDiagnostics.com/GramPRO)