



**Princeton BioMeditech Corporation**

*Makers of Advanced Rapid Diagnostic Products*

P.O. Box 7139, Princeton, NJ 08543-7139, USA

Tel: 732-274-1000

4242 U.S. Highway 1, Monmouth Junction, NJ 08852-1905, USA

Fax: 732-274-1010

[www.pbmc.com](http://www.pbmc.com) / [info@pbmc.com](mailto:info@pbmc.com)

---

September 15, 2021

MEMO

To: Customers of *Status* COVID-19/Flu test kit  
Re: Expiration dating

Dear Valued Customer,

The shelf life of the *Status*<sup>TM</sup> COVID-19/Flu A&B test kit ("Product") was initially established as 6 months due to the rapid and recent development of the Product and according to regulatory requirements. PBM has continued to conduct real-time stability testing to extend the expiration date. The Product in your possession may now have an expiration date that is later than the labeled date.

The affected Product lot numbers, originally labeled expiration date, and updated expiration date are listed below. These Products were originally labeled with a 6-month expiration date. All Products now have an 11-month expiration date. Please retain this memo with your Product records.

<b>Lot Number</b>	<b>Printed Expiration Date</b>	<b>11-month Expiration Date</b>
881B11	2021-08-31	2022-01-31
881B21	2021-08-31	2022-01-31

As additional data become available, we expect to communicate future improvements to you in additional memos.

We are pleased to communicate this improvement in the utility of the Product and look forward to assisting you further. Please contact us if you have any questions regarding this memo. E-mail: [techservice@pbmc.com](mailto:techservice@pbmc.com)

Sincerely,

---

Young-A Kim, Ph.D.  
Director, QA/RA  
PBM

---

Roger Kang  
President  
LifeSign, LLC